




# Sleeptracker-AI<sup>®</sup>

## Application User Guide






# Sleeptracker-AI<sup>®</sup>

Your **Tempur-Pedic<sup>®</sup> Sleeptracker-AI<sup>®</sup>** monitor and its companion application are powered by Sleeptracker technology to accurately and continuously monitor a sleeper's movements, heart rate, and respiration rate throughout the night.

Easy-to-understand sleep graphs visually distinguish between periods when you were awake, in light sleep, deep sleep, or REM sleep. These are quantified with a "Sleep Quality" metric. Goals and dynamic coaching encourage getting a better understanding of your sleep and ways you can improve your overall sleep health.




The Tempur-Pedic<sup>®</sup> Sleeptracker-AI<sup>®</sup> app lets you set a sleep cycle alarm that will help you Wake Up at the Optimal Time<sup>™</sup> in your sleep cycle so that you wake up refreshed and more energized throughout the day.

# Setting Up the Tempur-Pedic® Sleeptracker-AI® Monitor

To get started, make sure your Tempur-Pedic® Sleeptracker-AI® processor unit is powered on and the sensor cables are plugged in as indicated.


Next, install the Tempur-Pedic® Sleeptracker-AI® app onto your phone. The first time you use the app, it will guide you through the easy process of creating a user account and connecting your smart bed to your home Wi-Fi network.

Your smart bed must be connected to a Wi-Fi network in order to sync data with the Tempur-Pedic® Sleeptracker-AI® app.




# Home Page

The Home page shows an overview of daily progress toward your sleep goal. Important sleep statistics are viewable in the center ring. Daily tips, personal statistics, and comparisons from the Sleeptracker-AI® Sleep Coach are accessible at the bottom of the page — these tips will become more personalized as the app gets to know your sleep habits.



The center ring cycles through the time spent in each sleep phase:



Percentage of Sleep Goal: Set your initial sleep goal at a realistic time for you. Adjust from Menu>Sleep Goals & Schedule.

Select Analysis to view your sleep graph and sleep log.



### Home button:

Touch this button from anywhere else in the app to return to the Home page.

### Sleep Analysis Report:

View your full Sleep Report.

### Application menu:

Manage your alarm, sleep goal, user profile, and account settings. Add a second sleeper or get help using the app.

### Remote Control:

Control your adjustable base through the app.


### Air Quality:

Monitor air purity, CO2, humidity and the temperature of your bedroom.

# Sleep Recording

Tempur-Pedic® Sleeptracker® records your sleep automatically each night by default. Optionally, select the Moon button on the home page to start a recording when you are ready to fall asleep, versus laying in bed reading or watching TV, to improve the accuracy of calculated sleep efficiency values and for a more accurate Sleep Quality metric. Manual recording is also best for short naps you wish to record outside of your regular sleep routine.

Select the Moon button to begin a sleep recording and adjust your sleep cycle alarm settings.



The Status icon indicates your connectivity to the Sleeptracker-AI® cloud. Select the icon to view your connection status, add a second sleeper, change bed name, or change the side of the bed you sleep on.



# Sleep Cycle Alarm

The **Sleep Cycle Alarm** wakes you during your lightest phase of sleep in a given window of time. The alarm set time is the latest time at which the alarm will sound. Access the alarm settings from Menu > Sleep Cycle Alarm.

Set the latest time at which the alarm will sound.

**Sleep Cycle Alarm**

6:50AM

Repeat: Mon Tue Wed Thu Fri

S M T W T F S

Window

10min

The alarm will wake you at the optimal time between 6:40AM and 6:50AM.

**Alarm Type**

Vibrate Power Base

You'll use the the Sleeptracker® app on your phone or tablet to turn off the alarm. Note that the vibration may affect both sides of the bed.

Music on iPhone or iPad

Ocean Waves

HOME ANALYSIS MENU REMOTE AIR QUALITY

Choose an alarm window. For example: a 10 minute window will cause the alarm to fire between 6:40-6:50am when you are in a light phase of sleep.


Alarm Type: Select between a vibrating power base\*, music or a preset alarm sound to wake you up.

\* The Vibrate Power Base alarm type is only supported by compatible Tempur-Pedic® Power Bases



# Sleep Log

The sleep graph indicates periods of deep sleep, light sleep, REM sleep and awake time over the course of the night.



Analysis table of contents.

Tap on brief awake, sunrise or alarm icons for more info.

Heart rate graph

Respiration rate graph.


Scroll to the bottom of the page to Edit bed time or rise time.\*

\* For example: if you started your sleep recording too early, or if you manually start/stop your sleep recording and forgot to stop your recording in the morning





# Sleep Phases



**Awake:** fully conscious, alert and aware.

**REM:** characterized by rapid eye movement (REM) and increased cerebral blood flow relative to deep sleep. This is the cycle of sleep associated with dreaming.

**Deep sleep:** characterized by a lack of bodily movement and decreased cerebral blood flow, with associated changes in heart rate and respiration rate.

**Light sleep:** closer to wakefulness than deep sleep.